


Metabolic Monitor Sample Indirect Calorimetry Report

COSMED
Rome - ITALY
www.cosmed.com

Test date	Time
28/02/2018	12:30

Subject	John Doe	ID	--	Gender	Male	Age	51	D.O.B.	04/03/1967	Weight	102.0 Kg	Height	195.0 cm	BMI (kg/m2)	25.0
Predicted Set	For Comparison Harris-Benedict	Test Position	Undefined	Agitation	N/A	Sedation	N/A	Body Temp (°C)	0	Resting period	N/A	Fasting period	N/A	UN (g/day)	0

Indirect Calorimetry Report - Ventilator Test



REE		2634	RQ	0.87
125 %pred		Kcal/day	npRQ: --	--
V02	377	VCO2	327	Vt
mL/min	mL/min	mL/min	933	mL
Substrates			Rf	18.9
FAT	45	CHO	55	PRO
%	%	%	0	%
Variability			FiO2	40.03
V02	5	VCO2	7	AVG Interval
%	%	%	%	05:00
Operator:				Duration
Signature:				00:11:00
				min

Sample printout of a hypothetical Q-NRG+ patient REE results in a comprehensive format to facilitate metabolic assessment.

Abbreviations: REE – Resting Energy Expenditure; Predicted Set – Compared to a selected predictive equation; V0₂ – Oxygen Consumed; VCO₂ – Carbon Dioxide Produced; RQ – Respiratory Quotient (VCO₂/V0₂); Vt – Tidal Volume; Rf – Respiratory Frequency; FiO₂ – Fraction of Inspired Oxygen; CHO – Carbohydrates; PRO – Protein; AVG – Average; EE – Energy Expenditure, npRQ - non-protein respiratory quotient; mL/min - Milliliter/minute; min - minute

Q-NRG+ Technical Specifications

Product		
Part Number	C09092-02-99 (EMA), C09092-12-99 (North America), C09092-22-99 (A/NZ)	
Intended Use	Resting Energy Expenditure (REE) measurement on mechanically ventilated and spontaneously breathing subjects*	
Standard Packaging	Q-NRG+, USB cable, power cable, User Manual	
Test Kit (Single-use)	Flow-REE, FiO ₂ and FeO ₂ /CO ₂ sampling lines, FiO ₂ Vent Adapter, HME or standard filter	
Measurement Modes		
Ventilator	Standard	
Canopy Hood	Optional	
Face Mask	Optional	
Main Parameters	Range	Accuracy
VO ₂	10-1000 mL/min	±3% or 5mL/min
VCO ₂	10-1000 mL/min	±3% or 5mL/min
RQ	0-2.00	±5% or 0.04
REE	0-7200 kcal/day	±5% or 36 kcal/day
Flowmeter	Ventilator	Canopy/Mask
Type	Disposable Pneumotach (Flow-REE)	Bidirectional digital turbine
Flow Range	0.01 – 1.6 L/s	0.05 - 2 L/s
Accuracy	≤2% or 100mL/min @1-25 L/min	≤ 2% or 100mL/min @1-25 L/min
Resistance	2.3 cmH ₂ O s/L @ 1 L/s	<0.25 cmH ₂ O s/L @ 1 L/s
Calibration	Automatic via Internal Blower	With 3L calibration syringe (monthly)
Gas Sensors	O ₂	CO ₂
Gas Exchange Sampling	Micro Dynamic Mixing Chamber (patented)	
Type	Galvanic Fuel Cell (GFC)	Digital NDIR
Range	0-75%	0-10%
Accuracy	<0.05% Vol	<0.05% Vol
Resolution	0.01% Vol	0.01% Vol
Calibration	Automatic via gas cylinder (monthly)	

For more information on this medical device, please refer to User Manual.

Baxter Healthcare SA
Thurgauerstrasse 130
8152 Glattpark, (Opfikon), Switzerland
T +41 44 878 60 00

www.baxter.com

Report	
Export Modes	USB, Bluetooth®
Export Formats	PDF, CSV, XML
Hardware	
Display	10.1" Transmissive TFT LCD, 1024x600, 65k colors, capacitive touch screen
Power	Battery: Li-Ion "smart" (3 hours autonomy) Main: 100V-240V ±10%; 50/60Hz, 1.5A @100VAC, 1.0A @240VAC
Wireless Connectivity	Bluetooth (2.1 + EDR Class II - Range 10 m line-of-sight)
Wired Connectivity	1 USB-device (5kV galvanic-insulated), 2 x USB Host, RS-232, LAN
Weight & Dimensions	4.65 kg (10.3lb), 31x21x27cm (12.2x8.3x10.6in)
Environmental Ranges	Temp. +10°C to +35°C. Humidity: 5-93% (non condensing). Atmospheric pressure: up to 3011m
PC Software (optional)	
Languages	Italian, English, Spanish, French, German, Portuguese, Greek, Dutch, Turkish, Russian, Chinese (Traditional & Simplified), Korean, Romanian, Polish, Czech, Norwegian, Hebrew
OS Requirements	Windows 7, 8, 10
Security & Quality Standards	
MDD (93/42/EEC Class IIa), Safety (Class I IEC 60601-1), EMC (IEC 60601-1-2), Telemetry (ETSI EN 301 489-17)	

* This device is intended for the measurement of REE with some limitations in accordance with labeling, within the following population:

Ventilator: ventilated subjects > age 10 and 10Kg (22lb)

Canopy: spontaneously breathing subjects >15Kg (33lb)

Mask: spontaneously breathing subjects > age 6 and 10Kg (22lb)

Intended Use: This device is intended for the measurement of REE for spontaneously breathing and ventilated patients, with some limitations in accordance with labeling, within the following population: spontaneously breathing subjects >15 Kg (33 lb), when tested with the Canopy dilution technique, ventilated subjects > age 10 and 10 Kg (22lb), and spontaneously breathing subjects > age 6 and 10 Kg (22 lb), when tested with Face Mask. This device is not suitable for operating in presence of flammable anesthetic gases or gases other than O₂, CO₂, N₂ and water vapor. The device is to be used by physicians or by trained personnel under the responsibility of a physician. The device is not intended as a continuous monitoring device for surveillance of vital physiological processes. Warnings: This device measures clinical parameters used to aid diagnosis and it is intended only as an adjunct device in patient assessment. In case of disturbing conditions, the shutdown is allowed because the safety of the device towards patients and operators is not affected, since the final evaluation is performed on the outcome data measured during a complete test. No modification of this device is allowed. For more information on this medical device, please refer to User Manual.



Baxter is registered trademarks of Baxter International Inc.
Q-NRG+ is a registered trademark of COSMED.
GLBL/QNRG/19-0002 08/2019